

RAPID ANALYSIS OF CYSTINE FOR THE DIAGNOSIS OF RENAL CALCULI

Jean-Philippe VEYSSIER ¹, Tijen TANYALCIN ²

1-Biochrom, Science Park, Cambridge, UK, 2-Tanyalcin Medical Laboratory, Izmir, Turkey

AIM

This study presents a rapid screening method for cystine measurement in urine with high peak resolution and accuracy and improved run time; one hour per sample.

INTRODUCTION






Classical cystinuria is the most common inborn error of aminoacid transport. The fact that the cystine is the least soluble natural aminoacid, its over excretion frequently results to the formation of cystine calculi in renal pelves, ureters and bladder along with infection, obstruction and with high risk of recurrent urinary stone formation that even may direct to renal failure. Cystine stones are difficult to treat and requires lifelong therapy. Diagnosis by quantification of cystine in urine and monitoring the patient with cystine excretion have a great impact on this disorder.

METHOD


This short program we developed by using Biochrom 30 Aminoacid analyser enables 15 analysis to be performed a day. The separation is achieved using 20 X 4.6 mm physiological high resolution column using predominantly buffer CII (pH: 3.15). Norleucine is used as an internal standard. The program is shown below. Single aminoacid standards are run by the program to identify the peaks first, than master mix of a working standard as a mixture of amino acids, acidics and Neutrals A6407-sigma and Basics A6282-sigma with SSA-Norleucine is prepared and run by the program. The single peaks of single amino acid standards (500 µM) prepared in loading buffer and the working standard containing aminoacids 500 µM in final concentration and cystine 250 µM are shown below.

RESULTS

No.	Time	Temp	Buffer	Pump	Nin	Reac	Comments
1	15:00	40°C	25.0ml/h	ON	ON		
2	06:00	95°C	25.0ml/h	ON	ON		
3	06:00	95°C	25.0ml/h	ON	ON		
4	02:00	60°C	0	OFF	OFF		Reset
5	15:00	60°C	2	31.2ml/h	OFF	OFF	Load
6	06:00	95°C	25.0ml/h	ON	ON		
7	15:00	60°C	2	25.0ml/h	ON	ON	
8	06:00	95°C	25.0ml/h	ON	ON		
9	06:00	95°C	25.0ml/h	ON	ON		
10	02:00	60°C	0	OFF	OFF		Reset
11	15:00	60°C	2	31.2ml/h	OFF	OFF	Load
12	06:00	60°C	2	25.0ml/h	ON	ON	



Proline 440 nm



440 nm
L_proline 500 µM

METHOD COMPARISON



11 real patient urine samples with a different wide range of cystine results are assayed by both full amino acids program and short cystine program. The data was checked for distribution and it was non parametric.

	Tests of Normality			
	Kolmogorov-Smirnov	Shapiro-Wilk	Anderson-Darling	Prandis-1915
Urine_AA	.311	.11	.054	.727
Short_Cystine	.311	.11	.054	.727

* Lilliefors Significance Correction

Rank	N	Mean Rank	Sum of Ranks
Short_Cystine - Urine_AA	11	5.50	60.50
Urine_AA - Short_Cystine	11	6.00	66.00
Total	22		126.50

Test Statistics ^a	
W	Short_Cystine - Urine_AA
Z	Short_Cystine - Urine_AA
Asymp. Sig. (2-tailed)	1.000



The results of both programs showed a strong linear correlation as shown above by linear regression analysis

Bland-Altman Bias plot showed that all 11 points are between 1.96 SD

COMPARISON BETWEEN THE FULL AMINO ACID PROGRAM WITH SHORT CYSTINE PROGRAM.

The data has been considered by means of the mean median values obtained from participating labs in the ERNDUM Scheme 2010

FF testing material	Lab no 392	Cystine result in Full programme (µmol/L)	Scale Standard Deviations	Cystine result in Short Cystine Programme (µmol/L)
150	Number of labs enrolled = 149 Mean = 42.7 Median = 42.0 SD = 5.19	45.2	0.0 - 0.5 SD	52.43
151	Number of labs enrolled = 175 Mean = 63.8 Median = 67.0 SD = 8.15	63.32	(-0.5) - 0.0 SD	56.56
152	Number of labs enrolled = 173 Mean = 31.7 Median = 32.0 SD = 3.78	29.40	(-1.0) - (-0.5) SD	36.16

CONCLUSION

This short program will ease the work load of clinical laboratories that have a continuous requirement for the analysis of urine specimens from patients with renal stones to screen for cystinuria or follow up the patients under treatment. Concerning the cost effectivity, short cystine program will be a useful tool for effective time saving and efficient reagent usage.